Bleach 'n Smile Automix

Bleichgel

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Directions for Use

Bleach 'n Smile Automix is a bleaching gel intended solely for use in dental surgeries. It contains bleaching gel and activator in an auto-mixing 4:1 dual syringe (cartridge). One double syringe contains 5 g of materials. The ready-mixed gel contains 35% hydrogen peroxide.

Indications

- Brightening of non-vital teeth
- Brightening of discoloured teeth as an alternative to a crown or veneer
- Treatment of tooth discolorations caused by drugs (e. g. Tetracycline)
- Aesthetic tooth brightening

Contraindications

- Exposed dentine
- Defective filling margins
- Hypersensitive patients
- Allergy to or intolerance of the ingredients
- Pregnancy and lactation (due to insufficient data on the effects on pregnant and lactating women)

Many patients may experience slight sensitisation in their teeth during or shortly after treatment, but this normally subsides after a short time.

Stop treatment immediately if undue pain is being experienced.

Preparations

Determine the tooth colour by comparing it with the Vita® Shade Guide. To evaluate the colours, please re-sort the Vita® Shade Guide as follows

Values sorted by hue:

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Sorted according to brightness value:

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B2-B1 = Values are increased twofold

A3-A2 = Values are increased fourfold

Use air flow or a cleansing paste to clean the patient's teeth. Use a cheek/lip retractor or similar instrument. Protect the gingiva (e. g. with Bleach 'n Smile Dental Dam or by applying a rubber dam). Do not anaesthetise the patient.

The pure activator (red component) is caustic. The mixed bleaching gel is irritating to the eyes and mucous membranes. Wear gloves and goggles when using the material. Avoid contact with the eyes, skin and mucous membranes.

Application

Remove the cap and replace it with a mixing tip. Now squeeze the gel from the syringe until a homogeneously coloured mass with gel-like consistency is coming out of the tip. Discard the initial inhomogeneous mass.

Storage for long periods at high temperatures may sometimes cause bubbles to form in the gel. This can also lead to inhomogeneous and more liquid mixtures, especially in the last quarter of the dual syringe. If this occurs, do not use the affected parts of the material.

The gel should be applied in a thickness of about 1 - 2 mm to the labial surfaces of the teeth to be bleached. The point of the mixing tip is designed to spread the gel evenly.

To enhance the effect, you can use a standard plasma arc light or laser. **In this case, both the dentist and patient should wear safety goggles.** Please follow the manufacturer's instructions. Do not overheat.

After a reaction time of about 10 minutes, aspirate the gel and wipe off the residue. If the bleaching effect is insufficient, the procedure can be repeated twice until the desired effect has been achieved. Ensure beforehand that the gingiva has adequate protection. Then determine the tooth colour again.

If there is still enough material for another bleaching treatment, remove the mixing tip and replace the blue cap. Store the material in a refrigerator. Use a new mixing tip for the next treatment.

Note: Any composite fillings and restorative materials made of ceramic or metal are not bleached. Fillings should therefore be replaced 2 weeks after treatment.

Non-Vital Teeth

Extra gel is inserted via a lingual access. For this, the root canal filling is removed up to 2 mm beneath the cementenamel junction. Next, cover with Glasionomer cement while ensuring a correct marginal fit. Now carry out the brightening procedure as described above.

Finishing

After the final bleaching procedure, aspirate the gel, wipe it off and then rinse thoroughly with water. Now dry the tooth surface and polish with polishing paste. Next, remineralisation should be carried out with fluoride solution or gel. This is also suitable for treating pain if sensitisation occurs.

The patient should not eat, drink or smoke for 1 hour after treatment. Substances with strong colours (e. g. red wine, coffee) and foods containing acid should also be avoided in the days immediately following the treatment.

Shelf Life / Storage

Do not use Bleach 'n Smile Automix after the expiry date, because an adequate bleaching effect is no longer guaranteed.

Refrigerate the bleaching gel at 3 - 9 °C. Do not freeze.

The product contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and colouring.

Troubleshooting

Insufficient bleaching effect	Repeat the bleaching procedure (several times if necessary)
Patient is experiencing pain	The gel has come into contact with exposed dentine or a loose filling (please cover before treatment). The patient is particularly sensitive (depending on the individual tooth structure)
There are white spots on the enamel after the bleaching procedure	Local demineralisation or dehydration. Reversible. Disappears after fluoridation.
Gingiva is coloured white after bleaching (procedure can be reversed)	The gingiva is not adequately protected, or the protection is loose
The gel from the cartridge has an inhomogeneous colour (colourless/red)	The gel has air bubbles in it as a result of being too stored at too high a temperature or for too long
Reduced composite-enamel adhesion after bleaching procedure	Peroxide residue. Leave for 1 - 2 weeks before restoration.

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Schütz Dental GmbH • Dieselstr. 5-6 • 61191 Rosbach/Germany