BLEACH DAM

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BLEACH DAM

Light Cure Flowable Composite for Gingiva Protection

BLEACH DAM is a **visible light cure flowing microhybrid composite** for protecting the gingiva adjacent to teeth being whitened. It is based on methacrylate-resin and inorganic filler particles. The blue color facilitates optical control during application.

Indication for BLEACH DAM

- Gingiva protection during bleaching procedures

Application of BLEACH DAM

Isolate the gingiva with cotton rolls and gauze.

BLEACH DAM builds on the tooth surface a barrier to the gingiva.

Apply BLEACH DAM directly from the syringe using the enclosed needle tips at the gingival border on the enamel. The layer should have a thickness of approx. 1-1.5 mm.

Light cure BLEACH DAM with a commercially available light cure unit for **10 seconds** per arch.

After bleaching procedure is finished, remove BLEACH DAM with a spatula or tweezers.

Storage

Do not store above 25 °C (78 °F). Cool storage prolongs shelf life.

Additional Notes

- Do not use any resin to adjust viscosity of composite material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.